



Convert Gas Charged to...



Gas Charged Cartridge

Maintenance Free



Maintenance Free
 ECT Cellular Tube

Save time and money while reducing risk

- Reduced operating costs
- No more nitrogen charging
- Longer operating life
- Safer operations

Easy installation

- ECT retrofit kit preassembled
- Drop in place
- All-in-one easy installation
- Requires no vessel/piping modifications or special tools
- Lifting eye for easy removal and installation
- Retrofits in CoorsTek® products

Boost reliability and performance

- Ready to use upon installation
- No overcharging
- No premature failure
- No effect on pressure drop
- Operates at suction pressures to full vacuum

Performance Pulsation Control® is a registered trademark of Performance Pulsation Control, Inc. Coorstek® is a registered trademark of COORSTEK, INC. Performance Pulsation Control and its products are not affiliated with Coorstek. © 2019 Performance Pulsation Control, Inc. All rights reserved.

- **Convert to Maintenance Free** - PPC specifically developed the ECT cellular tube to convert gas charged to maintenance free
- **Stop Gas Charging** - ECT uses millions of nitrogen-filled, closed cell microcells to eliminate gas charging your product
- **Gain Consistent Performance** - ECT's rapid response and large surface area gives you consistent performance
- **Easy to Use** - ECT is an easy replacement - requires no piping/vessel modifications or special tools



Available ECT Retrofit Kits

Status Flow[®]

- SFVM Series
- SFT Series
- SFV Series
- SFM Series
- SFG Series
- SFU Series

CoorsTek[®]

- GM, MGM, SGM Series
- T, ST Series
- G, SG Series
- M, MM, SM Series

ECT System available in complete units.

Performance Pulsation Control[®] is a registered trademark of Performance Pulsation Control, Inc. Coorstek[®] is a registered trademark of COORSTEK, INC. Performance Pulsation Control and its products are not affiliated with Coorstek. © 2019 Performance Pulsation Control, Inc. All rights reserved.